

Pico mini Tattoo removal



Pico Mini Laser αφαίρεσης Tattoo.

Ένα ισχυρό compact Nd Yag Q-Switch Laser ιδανικό για αφαίρεση PMU-Tattoo.

ΚΩΔΙΚΟΣ: N/A | **Κατηγορίες:** [Laser - IPL](#), [Εξοπλισμός Tattoo - PMU](#) |

Φωτογραφίες



Περιγραφή Προϊόντος

Το Pico Mini είναι ένα Nd Yag Laser Q-Switch στα 1000 W και συχνότητα έως 15Hz ειδικό στην αφαίρεση ημιμόνιμου-μόνιμου τατουάζ.

Διαθέτει 3 κεφαλές:

- 1064nm : Αφαίρεση σκούρων χρωστικών
- 532nm : Αφαίρεση ανοιχτόχρωμων χρωστικών
- Carbon : Ανάπλαση

Compact desing για εύκολη μεταφορά και ανθεκτική, εργονομική κεφαλή αλουμινίου για άνεση κατά την εφαρμογή.

Laser Type	Suppise and ruby Q-switched laser instrument
Certification	Medical Ce IEC60601 TUV SPS
Wavelength	1064nm 532nm
Function	Pigment deposit dispersing, Tattoo removal, Skin care
Cooling System	Wind-cooler
Machine Size	35*45*157cm
Display	8" color touch screen
Pulse Energy	800mj
Width of Pulse	6ns
Frequency	1-6Hz
Instruction Light	Infrared ray indicator
Spot Diameter	1-8mm
Net Weight	24 kg
Voltage	220V/110V
Service	CEM/CEM Service
Language	English Standard, or Your Language to OEM
Delivery Time	3-7 Working Days
Warranty	1-2 Year, Technology Supports Forever

